

Product Datasheet

Eme1 antibody (orb76616)

Description	Mouse monoclonal antibody to Eme1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	Human EME1 (60-65Kd)
Target	EME1 (Msm4)
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	0.8mg/ml
Note	For research use only
Application notes	Sp2/0 Ag.14;
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q96AY2
Entrez	146956
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline use at 1:1000 for western blot
Expiration Date	12 months from date of receipt.