


Product Datasheet

EPCAM antibody (orb76621)

Description	Rat monoclonal antibody to EPCAM
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-Fr, IP, WB
Immunogen	Thymic epithelial cell line
Target	EpCAM
Preservatives	Affinity purified antiserum containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	X63Ag8.653;
Isotype	IgG2a kappa
Clonality	Monoclonal
Uniprot ID	Q99JW5
Entrez	17075
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at 0.5ug/ml for IHC on mouse material
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...