

Product Datasheet

FGFR4 antibody (orb76622)

Description	Mouse monoclonal antibody to FGFR4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Cytoplasmic domain of human FGFR4 fused to maltose binding protein
Target	FGFR4
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	NS1 myeloma
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O55111 , P22455
Entrez	116501
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try at WB: 1/100. Detects a band of approximately 105 kDa (predicted molecular weight: 88 kDa).
Expiration Date	12 months from date of receipt.

