



Product Datasheet

Ribosomal Protein L10L antibody (orb766231)


Description	Rabbit polyclonal antibody to Ribosomal Protein L10L
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RPL10L. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various cells
u...



Western
blot
analysis of
lysates from
Ju...



Western
blot
analysis of
the lysates
fro...