



Product Datasheet

Ribosomal Protein S4Y1 antibody (orb766250)


Description	Rabbit polyclonal antibody to Ribosomal Protein S4Y1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RPS4Y1. AA range:214-263
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	29
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western blot analysis of 293T mouse-sple...



Immunohistochemical analysis of paraffin...