

Product Datasheet

FXR2 antibody (orb76628)

Description	Mouse monoclonal antibody to FXR2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IP, WB
Immunogen	His6-FXR2s fusion protein produced in a bacterial pET system
Target	FXR2
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P51116
Entrez	23879, 9513
Dilution Range	Use at an assay dependent dilution. Detects a band of approximately 92 kDa.
Expiration Date	12 months from date of receipt.