


Product Datasheet

GAD65 antibody (orb76630)

Description	Mouse monoclonal antibody to GAD65
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Synthetic peptide: GSGFWSFGSEDGSGDSEN conjugated to KLH, corresponding to amino acids 4-22 of Human GAD65
Target	GAD65
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Sp2/0 Ag.14
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9UGI5
Entrez	14417 , 24380 , 2572
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of hum...