


Product Datasheet

Gemin3 antibody (orb76631)

Description	Mouse monoclonal antibody to Gemin3
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC, IP, WB
Immunogen	Recombinant 6His-tag COOH-terminal domain of Gemin3 (amino acids 368-548)
Target	Gemin 3
Preservatives	Purified antibody (from supernatant) containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0;
Isotype	IgG1
Clonality	Monoclonal
Entrez	11218
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try:WB: Use at an assay dependent dilution. Detects a band of approximately 105 kDa (predicted molecular weight: 99 kDa).Flow Cyt: Use 1µg for 106 cells.ICC/IF: 1/50.
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...