



Product Datasheet

Snk antibody (orb766347)


Description	Rabbit polyclonal antibody to Snk
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PLK2. AA range:291-340
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western Blot analysis of MCF7 cells usin...



Western blot analysis of lysates from He...