

Product Datasheet

HCG antibody (orb76640)

Description	Mouse monoclonal antibody to HCG
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	Purified human beta-HCG
Target	HCG beta
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	X63Ag8.653;
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P0DN86
Entrez	1082
Dilution Range	Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used in ELISA
Expiration Date	12 months from date of receipt.