



Product Datasheet

Synphilin-1 antibody (orb766411)


Description	Rabbit polyclonal antibody to Synphilin-1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Synphilin-1. AA range:797-846
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100
Expiration Date	12 months from date of receipt.



Western Blot analysis
of various cells u...



Immunohistochemical
analysis of paraffin...



Immunohistochemistry
analysis of Synphil...