


Product Datasheet

HDAC2 antibody (orb76642)

Description	Mouse monoclonal antibody to HDAC2
Species/Host	Mouse
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, FC, ICC, IHC-P, IP, WB
Immunogen	Synthetic peptide, corresponding to amino acids 473-488 of Human HDAC2
Target	HDAC2
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q92769 , P70288 , P56519
Entrez	15182 , 84577 , 395635 , 407223 , 3066
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1/1000 - 1/4000. Detects a band of approximately 55 kDa (predicted molecular weight: 55 kDa).
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell lysat...