


Product Datasheet

HDAC3 antibody (orb76643)

Description	Mouse monoclonal antibody to HDAC3
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Human C-terminal peptide (identical in mouse). Sequence: aa 411-428 (NEFYDGDHDNDKESDVEI)
Target	HDAC3
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	O15379, O88895
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline start at 1:400-1:800 for western blot. Detects a band of approximately 49 kDa (predicted molecular weight: 49 kDa).
Expiration Date	12 months from date of receipt.



Western blot analysis of NIH3T3 mouse fi...