

Product Datasheet

HEF1/ CasL antibody (orb76645)

Description	Mouse monoclonal antibody to HEF1/ CasL
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Fusion protein, corresponding to amino acids 82-398 of Human HEF1.
Target	HEF1/NEDD-9
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q14511
Entrez	4739
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;IP: Use a concentration of 1 µg/ml. Use at a concentration of 1 µg/200µg whole cell lysate.IF: 1:100
Expiration Date	12 months from date of receipt.