



Product Datasheet

TSLC1 antibody (orb766525)


Description	Rabbit polyclonal antibody to TSLC1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CADM1. AA range:393-442
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49
Expiration Date	12 months from date of receipt.



Western Blot analysis of various cells u...



Western blot analysis of 293T lysis usin...



Western blot analysis of 293T lysis usin...