


Product Datasheet

Histone H4 Me3K20 antibody (orb76654)

Description	Mouse monoclonal antibody to Histone H4 Me3K20
Species/Host	Mouse
Reactivity	A. thaliana, Drosophila, Human, Mammal, Other
Conjugation	Unconjugated
Tested Applications	ChIP, DOT, ELISA, WB
Immunogen	Synthetic peptide of Human Histone H4 (Immunogen sequence: 15-AKRHR(Kme3)VLRDN-25 (branched))
Target	Histone H4 (tri methyl K20)
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P84040, P59259, P02304, P62799, P62805
Entrez	3771938, 3772314, 8368, 8361, 3772319, 554313, 3772666, 8294, 3772254, 3772708, 8365, 8364, 8359, 3772129, 3772317, 3771941, 8370, 3772113, 3771893, 3772096, 41773, 447787, 3772509, 3771935, 8366, 318846, 121504, 8363, 3772325, 3772172, 3771908, 8362, 3772211, 3771947, 3772519, 8367, 3771854, 8360
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: Use at 1:4000. Predicted molecular weight: 11 kDa.
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela cell lysat...