


Product Datasheet

HIV1 P17 Gag antibody (orb76658)

Description	Mouse monoclonal antibody to HIV1 P17 Gag
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	Bacterially expressed, hexahistidine amino-terminal tagged HIV-1 p17 Gag protein (clade B, HXB-3 isolate).
Target	HIV1 p17
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2b kappa
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try: WB: 1/2000.
Expiration Date	12 months from date of receipt.



Western blot analysis of semi purified H...