



Product Datasheet

CHD1L antibody (orb766632)

Description	Rabbit polyclonal antibody to CHD1L
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CHD1L. AA range:545-594
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	101
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from N...



Western blot analysis of lysates from He...