



Product Datasheet

EF-1beta antibody (orb766634)

Description	Rabbit polyclonal antibody to EF-1beta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human EF-1beta. AA range:83-132
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from K...



Western blot analysis of lysates from K5...