



Product Datasheet

SAP 155 antibody (orb766643)


Description	Rabbit polyclonal antibody to SAP 155
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human SAP 155. AA range:75-124
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	145
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from K...



Western blot analysis of lysates from Ju...



Immunohistochemical analysis of paraffin...