



Product Datasheet

A1BG antibody (orb766662)


Description	Rabbit polyclonal antibody to A1BG
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from A1BG . at AA range: 100-180
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	28
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from K...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...