



Product Datasheet

TRPV3 antibody (orb766700)


Description	Rabbit polyclonal antibody to TRPV3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the Internal region of human TRPV3.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90
Expiration Date	12 months from date of receipt.



Western Blot analysis of extracts from C...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...