



Product Datasheet

RRP7A antibody (orb766718)

Description	Rabbit polyclonal antibody to RRP7A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human RRP7A.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	32
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562, 293, SK-N-BE(2)C, and mouse eye.



Immunohistochemical analysis of paraffin-embedded tissue.