


Product Datasheet

SCP-2 antibody (orb766721)

Description	Rabbit polyclonal antibody to SCP-2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the Internal region of human SCP-2.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	175
Expiration Date	12 months from date of receipt.



Western Blot analysis of SKOV3, BT474, J...



Immunohistochemical analysis of paraffin...