



Product Datasheet

ACTA2 antibody (orb766725)


Description	Rabbit polyclonal antibody to ACTA2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human α -SMA.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42
Expiration Date	12 months from date of receipt.



Western Blot analysis of MFC cells using...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...