



Product Datasheet

DR3 antibody (orb766740)


Description	Rabbit polyclonal antibody to DR3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from DR3 . at AA range: 230-310
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using DR3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using DR3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using DR3 antibody.