



Product Datasheet

IL10 antibody (orb766743)

Description	Rabbit polyclonal antibody to IL10
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL10. AA range:71-120
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20
Expiration Date	12 months from date of receipt.




Western Blot analysis of HepG2, Mouse li...



Immunofluorescence analysis of A549. 1,

...



Immunohistochemical analysis of paraffin...