



Product Datasheet

FAF1 antibody (orb766754)


Description	Rabbit polyclonal antibody to FAF1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human FAF1. AA range:461-510
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using FAF1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using FAF1 antibody.



Western blot analysis of lysate from K562 cells using FAF1 antibody.