



Product Datasheet

CD66e antibody (orb766777)


Description	Rabbit polyclonal antibody to CD66e
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CEACAM5. AA range:481-530
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76
Expiration Date	12 months from date of receipt.



Western Blot analysis of MCF7 cells using CD66e antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CD66e antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CD66e antibody.