



Product Datasheet

FGF-9 antibody (orb766784)


Description	Rabbit polyclonal antibody to FGF-9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FGF9. AA range:141-190
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23
Expiration Date	12 months from date of receipt.



Western Blot analysis of SW480 cells using FGF-9 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using FGF-9 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using FGF-9 antibody.