

Product Datasheet

HSV-1 envelope glycoprotein D antibody (orb76679)

Description	Mouse monoclonal antibody to HSV-1 envelope glycoprotein D
Species/Host	Mouse
Reactivity	HSV-1
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	HSV-1
Target	HSV-1 envelope glycoprotein D
Preservatives	Purified antibody containing PBS + 10mM sodium azide
Concentration	1.65mg/ml
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q69091
Entrez	2703444
Dilution Range	Optimal dilutions should be determined by the end user
Expiration Date	12 months from date of receipt.