

Product Datasheet

Human Milk Fat Globule Protein antibody (orb76680)

Description	Mouse monoclonal antibody to Human Milk Fat Globule Protein
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P, WB
Immunogen	Human milk fat globule membrane tissue / cell preparation
Target	Human Milk Fat Globule Protein
Preservatives	Purified antibody in PBS containing 10mM sodium azide
Concentration	2mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q08431
Entrez	4240
Dilution Range	Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used
Expiration Date	12 months from date of receipt.