



Product Datasheet

P2X3 antibody (orb766816)


Description	Rabbit polyclonal antibody to P2X3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human P2RX3. AA range:81-130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using P2X3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using P2X3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using P2X3 antibody.