



Product Datasheet

FOXP1 antibody (orb766818)


Description	Rabbit polyclonal antibody to FOXP1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FOXP1. AA range:628-677
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using FOXP1 antibody.



Western Blot analysis of customer's samples using FOXP1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using FOXP1 antibody.