



Product Datasheet

DPYD antibody (orb766832)


Description	Rabbit polyclonal antibody to DPYD
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human DPYD. AA range:351-400
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using DPYD antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using DPYD antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using DPYD antibody.