



Product Datasheet

SNAP 23 antibody (orb766836)


Description	Rabbit polyclonal antibody to SNAP 23
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human SNAP23. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using SNAP 23 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using SNAP 23 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using SNAP 23 antibody.