



Product Datasheet

SPTLC1 antibody (orb766837)


Description	Rabbit polyclonal antibody to SPTLC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from SPTLC1 . at AA range: 411-460
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using SPTLC1 antibody.



Western blot analysis of 293T SH-SY5Y cells using SPTLC1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using SPTLC1 antibody.