



Product Datasheet

Integrin alpha 4 antibody (orb766853)


Description	Rabbit polyclonal antibody to Integrin alpha 4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human ITGA4. AA range:571-620
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using Integrin alpha 4 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Integrin alpha 4 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Integrin alpha 4 antibody.