


Product Datasheet

Kallikrein 2 antibody (orb76686)

Description	Mouse monoclonal antibody to Kallikrein 2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-P
Immunogen	Balb/c tumor cell line transfected with human K2
Target	Kallikrein 2
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Spleen cells from immunized BALB/c.PSA transgenic mice were fused with the P3 myeloma cell line
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P20151
Entrez	3817
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;IHC-P: Use at an assay dependent dilution. As a guideline use at 18ug/ml for sections. Perform heat mediated antigen retrieval via the microwave method before commencing with IHC staining protocol.
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...