


Product Datasheet

Karyopherin alpha2 antibody (orb76687)

Description	Rat monoclonal antibody to Karyopherin alpha2
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Synthetic peptide, corresponding to amino acids 116/515 of Human KPNA2
Target	KPNA2
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	P3-x63Ag.8.653
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P52293 , P52292
Entrez	85245 , 16647 , 3838
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline start at 1:1000 for westerns
Expiration Date	12 months from date of receipt.



Western blot analysis of BL41