



Product Datasheet

CD201 antibody (orb766892)


Description	Rabbit polyclonal antibody to CD201
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from Endothelial protein C receptor at AA range: 141-190
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using CD201 antibody.



Western blot analysis of Rat-barin lysate using CD201 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CD201 antibody.