



Product Datasheet

GDI-2 antibody (orb766922)


Description	Rabbit polyclonal antibody to GDI-2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human GDI-2.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549, HT29, K562...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...