



Product Datasheet

HPA1 antibody (orb766958)


Description	Rabbit polyclonal antibody to HPA1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human HPSE. AA range:241-290
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	62
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using HPA1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using HPA1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using HPA1 antibody.