


Product Datasheet

Lamin A/C antibody (orb76696)

Description	Mouse monoclonal antibody to Lamin A/C
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant bacterially expressed full length protein
Target	Lamin A + C antibody [JoL5] - Nuclear Envelope Marker
Preservatives	Tissue culture supernatant containing 0.1% sodium azide
Note	For research use only
Isotype	IgG1 Kappa
Clonality	Monoclonal
Uniprot ID	P48678, P02545
Entrez	16905, 4000
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; ICC/IF: 1/10WB: 1/50. Detects a band of approximately 70 kDa (Lamin A) and 60kDa (Lamin C).
Expiration Date	12 months from date of receipt.



Western blot analysis of JoL5 cell lysat...