


Product Datasheet

Lamin A/C mutant R482W antibody (orb76699)

Description	Mouse monoclonal antibody to Lamin A/C mutant R482W
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	Lamin A Ig-fold; human (Immunogen sequence: 430-545 (6xHis-tagged fusion protein))
Target	Lamin A+C mutant R482W
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63-Ag8.653
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: 1:32,000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cell lysates...