


Product Datasheet

VDAC1 antibody (orb766995)

Description	Rabbit polyclonal antibody to VDAC1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human VDAC1. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31
Expiration Date	12 months from date of receipt.



Western blot analysis of lysate from K562...