



Product Datasheet

IL8R beta antibody (orb767002)


Description	Rabbit polyclonal antibody to IL8R beta
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human CXCR2. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEc cells using IL8R beta antibody (orb767002).



Immunohistochemical analysis of paraffin-embedded tissue sections using IL8R beta antibody (orb767002).



Immunohistochemical analysis of paraffin-embedded tissue sections using IL8R beta antibody (orb767002).