



Product Datasheet

CD169_Sialoadhesin antibody (orb767004)

Description	Rabbit polyclonal antibody to CD169_Sialoadhesin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human SIGLEC1. AA range:1321-1370
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	190
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells usin...



Western blot analysis of lysate from K56...