



Product Datasheet

CD163b antibody (orb767011)


Description	Rabbit polyclonal antibody to CD163b
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CD163L1. AA range:391-440
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using CD163b antibody.



Immunohistochemical analysis of paraffin-embedded tissue using CD163b antibody.



Immunohistochemical analysis of paraffin-embedded tissue using CD163b antibody.