



Product Datasheet

BMAL1 antibody (orb767028)


Description	Rabbit polyclonal antibody to BMAL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human BMAL1 around the non-acetylation site of Lys538. AA range:501-550
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	69
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...